

Einstein™ Site Prep Guide



Table of Contents

Document Information	3
Legal Notice.....	3
History of Changes	4
How to Use This Guide	4
User Information	5
Product Overview	5
Customer Responsibilities	5
Einstein Process Flow-Recommendation	6
Site Information	7
Environmental Conditions Einstein	7
Air Quality Considerations	7
Einstein Footprint	7
Product Safety	8
Einstein Printer	9
Printer Technical Specifications	10
Software	11
Overview Envision One RP Software®	11
System Requirements Envision One RP	11
Install Envision One RP Software	11
License Envision One RP Software	12
Accessories	13
Bottle Roller	13
Washing Unit.....	14
Curing Unit	15
Shopping List	17
Primary Supplies	17
Secondary Supplies	18
Appendix A – Site Prep Checklist	19
General Items	19
Checklist Signoff.....	20
Manufacturer	21

Document Information

Legal Notice

This document contains information that is confidential and proprietary to Desktop Metal, Inc. and/or its affiliates, including without limitation EnvisionTEC US LLC and ExOne Operating, LLC (each a “DM Company,” and together the “DM Companies”). This information is provided only to authorized representatives or customers of the DM Companies, and solely for the purpose of facilitating the use of DM Companies’ products and services. This document and its contents shall not be used or distributed for any other purposes or communicated, disclosed or copied except as agreed by a DM Company in writing. The information contained herein is provided for reference only and subject to change without notice. This document provides general information about the products described herein and is not a substitute for the Instructions For Use and Safety Data Sheets for said products. The DM Companies shall not be liable for omissions or for technical or editorial errors contained herein or for any damages whatsoever arising in connection with the furnishing or use of this document. This information is not intended to be used to determine the suitability or reliability of the user’s specific applications or environments; these determinations are the sole responsibility of the user, and the DM Companies disclaim all liability associated therewith. Without limiting the foregoing, the user is solely responsible for the use and operation of the products and services, including the disposal of waste products in connection therewith. This document does not supplement, replace or otherwise modify the terms and conditions that govern the purchase and sale or use of DM Companies’ products or services. Furthermore, nothing herein shall constitute a warranty; the only warranties for DM Companies’ products and services are those set forth in the express warranty statement in the terms and conditions of sale for said products and services. Desktop Metal, the DM Logo, Bound Metal Deposition, BMD, Live Parts, Studio System, Shop System, Fabricate, Fiber, Production System, Desktop Health, Desktop University, Flexcera, Einstein and ETEC are trademarks of Desktop Metal, Inc. EnvisionTEC, Envision One, cDLM, Vida, Perfactory, D4K, Hyperprint and Xtreme 8K are trademarks of EnvisionTec GmbH and its affiliates. ExOne, the ExOne Logo, Innovent, Innovent+, InnoventX, X25Pro, X160Pro, X1, S-Max, S-Print, CleanFuse, NanoFuse, and HydroFuse are trademarks of ExOne Operating, LLC or its affiliates. All other trademarks used herein are the property of their respective owners.

© 2022 Desktop Metal, Inc. All rights reserved.

History of Changes

Date	Changes	Version
Feb-2022	Document creation	1.0
March-2022	<ul style="list-style-type: none"> Added Legal Notice Updated Document Footer Updated Product Overview Updated How to Use This Guide Updated Einstein Footprint 	1.1
August- 2022	<ul style="list-style-type: none"> Added History of Changes Added Document Information Added User Information Added Site Information Updated How to Use This Guide Updated Environmental Conditions Einstein Updated Einstein Footprint Updated Printer Technical Specifications Updated Bottle Roller Updated Washing Unit Updated Curing Unit Added Manufacturer 	1.2

How to Use This Guide

This document serves as a comprehensive guide to prepare your site for the Einstein 3D printer.

User Information

Product Overview

Meet the Einstein 3D printer series, brought to you by Desktop Health. Meticulously designed for dental professionals, Einstein hits on all the key features essential to superior 3D printing: accuracy, speed, and versatility.

The secret to Einstein is in its DNA. Proprietary NanoFit 385 technology ensures that every application features stunning clarity, perfect fit, and an impressively natural-looking finish.

Equipped with Hyperprint™ technology, harnessing the power of heat and a closed-loop software upgrade, the Einstein series allows you to fabricate a variety of dental applications, from models to dentures and everything in between, with ease and at up to 50% faster speed than its predecessor.

Transform your workflow from fit to finish with Einstein. Available as Einstein, Einstein Pro, and Einstein Pro XL. Visit health.desktopmetal.com for more information.

Customer Responsibilities

NOTICE

Risk of equipment damage, risk of injury, and void warranty: comply with the below responsibilities.

Prior to installation, the customer is required to ensure the site is ready for the Einstein. These responsibilities are:

1. Provide a location to receive crated packages per individually listed shipping dimensions.
2. Provide a clear pathway from the receiving site to the installation site.
3. Provide resources to move the equipment from the receiving site to the installation location, with respect to:
 - a. The dimensions and weight of each system component.
 - b. Thresholds and corners.
4. A mover is required for unloading, uncrating, and moving the printer into location.
 - a. The printer requires a forklift or pallet jack to move the pallet to the site.
 - b. Carefully disassemble the crate.
 - c. Lift the printer out of the crate and place it in its final destination.
5. Create a component layout based on the system requirements included in this document.
6. Select a location that provides all services listed on each component page as necessary for operation.
 - a. Distance the printer from climate control systems. Do not direct airflow onto or within the printer's footprint. This includes air conditioning (window, wall, or similar), fans (ceiling, tabletop, or similar), heaters (space heaters, or similar).
 - b. Do not expose the printer to any direct sunlight or UV radiation
 - c. Keep the printer and accessories away from water sources such as sinks, tubs, wash stations, etc.
 - d. Do not place the printer near any vibration source, or on a wooden surface that can easily transfer vibrations.
 - e. Do not place the printer near polishing stations.
7. Meet all National/Governmental or Environmental Health and Safety (EHS) regulations.
 - a. See [Product Safety](#).
 - b. Review material storage as well as processing/operation guidelines.
8. Install all safety equipment required for operation, according to codes and regulations.
9. Perform risk analysis and address deficiencies.
10. Operate the Einstein indoors, as designed.

11. Locate the Einstein away from equipment that produces dust or vibration.



Note: If you have questions, please contact [Technical Support](#).

Einstein Process Flow-Recommendation

Process Steps

Model Processing:

1. Prepare the .stl file in CAD software.
2. Import the .stl file into [Envision One RP Software](#).
3. Use automatic placement to place parts.
4. Generate supports, if required.
5. Generate job in Envision One RP.
6. Send the build job to the printer.

At the Printer:

1. Install the clean material tray in the printer.
2. Prepare the material to print.
3. Add material to the material tray.
4. Install the clean build platform in the printer.
5. Close the printer hood.
6. On the touchscreen, navigate to the required build job.
7. Start the job. Wait for the job to finish.
8. Remove the build platform with printed model(s) from the printer.
9. Use the paint scraper to remove the model(s) from the build platform.

Cleaning Station:

1. Remove supports, if required.
2. For medical materials, follow instructions in the [Instructions for Use \(IFU\)](#) to clean models.
3. For non-medical materials, follow instructions in the [Material Best Practice guide](#) to clean models.

Post Curing:

1. Bring the clean and dry part(s) to the curing unit and place it in the unit.
2. Set the timer as per the recommended cure times (material dependent).

Finishing:

1. If part(s) were printed on supports, then sand away all support remains starting from a coarser grit to a finer one until the supported surface is smooth.
2. Final finishing processes are material dependent.



Note: The Einstein printer, curing unit, and washing unit are provided by Desktop Health. All other items to be purchased from third party suppliers. See [Shopping List](#) for recommendations.

Site Information

Environmental Conditions Einstein

Distance the printer from climate control systems. Do not direct airflow onto or within the printer's footprint. This includes air conditioning (window, wall, or otherwise), fans (ceiling, tabletop, or similar), heaters (space heaters, or otherwise).

- Do not expose the printer to direct sunlight or UV radiation.
- Do not place the printer near a doorway.
- Avoid temperature changes.
- Keep the printer and concentrator away from water sources such as sinks, tubs, wash stations, etc.
- Do not place the printer near polishing stations.

Ensure the room temperature remains within the following range:

- Minimum temperature of 23°C (73°F)
- Maximum temperature of 28°C (85°F)

Follow the temperature instructions for each material.

- These instructions can be found on the label of the material bottle.
- Extended instructions are located in the [Instructions for Use \(IFU\)](#) for medical materials, and the [Material Best Practice guide](#) for non-medical materials.



Important: The humidity should be below 45%.

Air Quality Considerations

Use the Einstein in a clean atmosphere. High levels of impurities in the air may find their way into printed parts.

- The room must allow heat generated from the system to dissipate at 0.03 m³/min airflow.
- Air conditioning units should be at least 3 m (10 ft) away from the printer with no airflow pointing directly at the printer.

Maintain a space with little or limited:

- Dust
- Smoke
- Steam

Einstein Footprint

NOTICE

Incorrect positioning of the printer can cause damage to the printer and/or to the printed models. Please follow the positioning guidelines.

Printer: 36.3 kg (80 lb).

Dimensions: 40.64 x 40.64 x 69.85 cm (16 x 16 x 27.5 in).

The Einstein must be installed on a sturdy surface that can comfortably hold over 36.3 kg (80 lb).

The printer requires a minimum amount of space for operation, maintenance, and heat dissipation.

- A minimum of 46 cm (18 in) behind the printer.
- A minimum of 61 cm (24 in) on the left and right side of the printer.
- A minimum of 122 cm (48 in) above the table-top.

Fig. 1 Einstein Clearance Front View





Fig. 2 Einstein Clearance Right Side View



Product Safety

Desktop Health's Einstein has been designed to comply with safety standards.

See the [Einstein Compliance and Safety Guide](#) > Product Safety Warnings for a description of all safety warnings associated with the Einstein printer.

General Safety

All general workplace safety rules should be followed when operating the printer. It is the responsibility of the user to ensure compliance with all local, regional, and national regulations. Additionally, it is the responsibility of the user to ensure that the system is installed and maintained properly by Desktop Health.

Material Safety

Safety data sheets (SDS) for materials used in the printing process are available in the Desktop Health Knowledge Base: [Safety Data Sheets](#). Read and understand the information provided in these documents prior to attempting to operate the printer or handle any media.

Einstein Printer

Fig. 3 Einstein Front View




Fig. 4 Einstein Back View





Fig. 5 Einstein Side View, Open Hood



Printer Technical Specifications

Unit Dimensions	Hood closed: 40.64 x 40.64 x 69.85 cm (16 x 16 x 27.5 in) Hood open: 40.64 x 47.76 x 84.77 cm (16 x 18.80 x 33.38 in)
Shipping Dimensions (estimate)	58.42 x 88.9 x 109.22 cm (23 x 35 x 43 in)
Unit Weight	36.3 kg (80 lb)
Shipment Weight (estimate)	74.84 kg (165 lb)
Electrical	Printer: 100-240V, 50/60Hz, 1 ph, 5A max
Ethernet (Not Required)	Cable: Cat 6 or better
USB	Type: USB 2.0 Standard A

Software

Overview Envision One RP Software®

Envision One RP automatically orients your model, adds supports, if necessary, and sends the file to the printer, resulting in your three-dimensional model. Everything that is printed using Desktop Health printers must pass through this software successfully.

System Requirements Envision One RP



Note: Envision One RP requires Windows Operating System.
The software is not compatible with macOS.

Operating System	Windows 10 or higher
Working Memory	>= 8GB RAM
Hard Drive	400 MB Free space
CPU	Multi Core Processor e.g. Core i5, >= 3GHz, >= 6MB Cache
Graphics	Dedicated 3D graphics card with >= 1 GB memory and OpenGL 4.3 and higher. It is recommended to use the Nvidia graphics card.

Install Envision One RP Software

It is recommended to install Envision One RP software prior to the arrival of your Einstein printer.

1. [Download Envision One RP Software](#).
2. Open the **Downloads** folder on your computer.
3. Find a file named *EnvisionOneRP.exe* and double-click it to start the installation.
→ A setup window opens.
4. Follow instructions on the screen to complete the installation.

License Envision One RP Software

Licensing Envision One RP provides access to the software's full functionality. There are two license types available:

- The **Standard license** allows you to use the whole software functionality for 12 months. After its expiration, you will have to request a new license.
- The **Trial license** allows you to use the whole software functionality for 14 days. After its expiration, you will have to request the **Standard license**.



Important: It is not possible to use the software without a license. If no license is available, the software is blocked.

See the Desktop Health Knowledge Base for further [licensing instructions](#), [software knowledge articles](#), and the [Envision One RP User Guide](#).


Accessories

Bottle Roller

Dual Motion Bottle Roller

The Dual Motion Bottle Roller is the recommended material mixing solution for all Desktop Health materials.

Fig. 6 Dual Motion Bottle Roller




Unit Dimensions	47 x 26 x 12 cm (18.5 x 10.24 x 4.72 in)
Shipping Dimensions	63.5 x 50.8 x 50.8 cm (25 x 20 x 20 in)
Unit Weight	5.1 kg (11.25 lb)
Shipment Weight	6.8 kg (15 lb)
Electrical Requirements	100-240V AC, 50Hz/60Hz

Washing Unit

PWA 2000

The PWA 2000 Parts Washing Apparatus is a brushless, contactless agitator and is the recommended washing unit for non-medical materials.

Fig. 7 PWA 2000, Open




Unit Dimensions	32 x 33 x 24.6 cm (12.60 x 12.99 x 9.69 in)
Shipping Dimensions	40.64 x 40.64 x 40.64 cm (16 x 16 x 16 in)
Unit Weight	4.6 kg (10.2 lb)
Shipment Weight	9.07 kg (20 lb)
Electrical Requirements	100-230V AC, 0.3A, 50Hz/60Hz External power adapter: 12V AC, 2A

Curing Unit

Otoflash

The Otoflash is the recommended curing unit for medical materials.


Fig. 8 Otoflash, Open



Unit Dimensions	31 x 31 x 14 cm (12.2 x 12.2 x 5.5 in)
Shipping Dimensions	38.1 x 33.02 x 17.78 cm (15 x 13 x 7 in)
Unit Weight	6.3 kg (13.9 lb)
Shipment Weight	8.16 kg (18 lb)
Electrical Requirements	100-230V AC, 15A, 50Hz/60Hz, 275 Watts

PCA 4000

The PCA 4000 Parts Curing Apparatus is the recommended curing unit non-medical printing materials.

Fig. 9 PCA 4000

Unit Dimensions	36.4 x 40.4 x 31.9 cm (14.33 x 15.9 x 12.55 in)
Shipping Dimensions	20 x 20 x 20 in. (50.8 x 50.8 x 50.8 cm)
Unit Weight	8.03 kg (17.7 lb)
Shipment Weight	11.34 kg (25 lb)
Electrical Requirements	100-230V AC, 4-6A, 50Hz/60Hz, 300 Watts

Shopping List

Primary Supplies

Primary supplies should be acquired prior to printer delivery. Obtain the proper tools to prepare for successful printing and post-processing.

- **99% Isopropyl alcohol (IPA):** breaks down uncured material from printed models. This product can be purchased from many big box stores, and from medical and business supply stores. Do not use denatured alcohol. Do not use less than 99%.
- **Absorbent paper towel:** used to clean the printer, material tray, and build platform. Inter-folded paper towel, found in labs or offices, is not as effective as a standard roll of absorbent paper towel.
- **Air compressor/forced air system:** required for drying models during post processing. Forced air also removes the 99% isopropyl alcohol from the surface of the models. Do not use canned air. For those without a forced air system, compressors can be purchased and installed into the post-processing station.
- **Desktop Health's Dual Motion Bottle Roller:** provides a gentle but highly efficient rocking and rolling motion, for 3D printing materials. The Dual Motion Bottle Roller almost eliminates the separation of resin and allows for an even distribution of the material. [Contact your Desktop Health Sales Representative](#) for more information.
- **Cone-shaped paint filters (included in the Starter Kit):** used to strain material. The filter removes particles that may be too small to see but may affect print quality. They will be used while pouring material from the tray into the material storage containers for maintenance. Each filter is approved for one-time use only. Do not leave filters in the bottle of material for longer than 5 minutes. Do not use coffee filters.
- **Desktop Health's curing unit:** 3D printed models must be post-cured to achieve the final, end-use state. Desktop Health offers the Otoflash curing unit for medical materials, and the PCA 4000 curing unit for non-medical materials. See Desktop Health's Knowledge Base for more information on the [Otoflash](#) and the [PCA 4000](#).



Note: Desktop Health only supports Desktop Health curing units.

- **Desktop Health's washing unit:** 3D printed models must be washed prior to post-curing. Desktop Health offers the PWA 2000 washing unit for non-medical materials. See Desktop Health's Knowledge Base for more information on the [PWA 2000](#).
- **Digital calipers (included in the Starter Kit):** required during printer calibration. The printer may need to be calibrated when new material is purchased, or for maintenance reasons. Digital calipers are used to check the accuracy of printed models against the computer's measurements. Find calipers from major manufacturers only. Look for accuracy of ± 0.02 mm.
- **Ethernet cable:** connects the printer to the local network. The cable inserts into the I/O panel behind the printer. The other end plugs into the network or operating computer. The printer can be operated by other Windows computers on the network if the CAT 6 is plugged into a modem or into a network outlet. If it is directly connected to the operating computer, it can only be accessed by the connected printer.
- **Garbage can:** for cleanup. Having a can on hand speeds up throughput, increasing efficiency. This needs to be easily accessible.
- **Nitrile gloves:** required when working with the printer or prior to touching anything that may come into contact with uncured material. See the [Safety Data Sheet](#) for proper handling guidelines.
- **Rubber spatula (included in the Starter Kit):** used to gently mix material in the material tray.
- **Storage containers for material:** used to store material after it is removed from the material tray. To avoid contamination, never pour used material back into the original bottle. Storage containers

must be completely clean, opaque, and have a strong seal to keep material fresh. Containers must not let any light through. See the [Safety Data Sheet](#) for proper handling guidelines.

- **USB drive (included in the Starter Kit):** a back-up for transferring print information from the operating computer to the printer. Between 8 and 32 GB size is recommended.

Secondary Supplies

Secondary supplies are recommended for the 3D printing workflow and may be acquired after printer delivery.

- **Magnification loupe:** used to look at the fine details of printed models.
- **Makeup or quality paint brush:** used to clean model(s) during post processing. Do not use a toothbrush as they are too rigid and can scratch printed models.
- **Metric Allen key set:** used to perform calibrations as directed by a Desktop Health certified technician.
- **Phillips screwdriver:** used to open the printer's front and back covers. The cover should only be opened when instructed by an authorized technician.
- **Post-processing tools:** Electric hand-tools, files, rasps, needle-nose pliers, sandpaper (220 grit +), and snips are all helpful tools for removing supports and polishing supported surfaces.
- **UV safety glasses (included in the Starter Kit):** used during printer calibration. They are not required for day-to-day operation.



Note: [Contact your Desktop Health Sales Representative](#) for more information on any Shopping List item.

Appendix A – Site Prep Checklist

Contact Information:	
Company name:	
Street address:	
City:	
Postal Code:	
Country:	
Phone:	
Printer Operator:	
Email:	
Phone:	
IT administrator:	
Email:	
Phone:	

General Items

Receiving & Site Access

- Receiving location meets requirement for size and weight of equipment ----- Y / N
- Clear path to install location (clear corners, doors, etc.) ----- Y / N
- Fork lift or pallet jack available to move items ----- Y / N
- Two individuals to lift 36.3 kg (80 lb) printer ----- Y / N
- Room accessible and Desktop Health staff access permitted ----- Y / N

General Facility Requirements

Minimum footprint

- Space available meets minimum footprint requirements ----- Y / N
- Room meets ventilation requirements for safe operation ----- Y / N
- Room is able to maintain 23-28°C (73-85°F) operational temperature ----- Y / N
- No vibration that will impact performance ----- Y / N
- Room is able to maintain less than 45% humidity ----- Y / N
- All necessary corporate, local, regional, or national approvals for operation ----- Y / N

Room Layout

- Required component clearances available 61 cm (24 in) on right and left sides ----- Y / N
- Required component clearances available 122 cm (48 in) above the tabletop ----- Y / N
- Required component clearances available 46 cm (18 in) behind the printer ----- Y / N

Safety

Completed a preinstall safety review ----- Y / N

Installed fire extinguishers s required by local code ----- Y / N

Printer

Power 100-240V, 50/60Hz, 1 ph, 5A max ----- Y / N

Materials

Printing material arrived on site and is undamaged ----- Y / N

Instructions for Use or Material Best Practice located and understood, mixing procedure understood, SDS located and understood ----- Y / N

Checklist Signoff

Sign below to acknowledge that the checklist is complete, accurate, and the site is ready for install:

Name of signatory: _____

Signature: _____

Title: _____

Date: _____



Note: If you have questions, please reach out to [Technical Support](#).

Manufacturer

Desktop Metal, Inc.
63 3rd Avenue
Burlington, MA 01803
www.desktopmetal.com

Desktop Health
c/o Desktop Metal, Inc.
63 3rd Avenue
Burlington, MA 01803
www.health.desktopmetal.com

EnvisionTEC US LLC (ETEC)
15162 Commerce Dr. S
Dearborn, MI 48120
www.etec.desktopmetal.com

EnvisionTec GmbH
Brusseler Str. 51
45968 Gladbeck
Germany

ExOne Operating, LLC
127 Industry Boulevard
North Huntingdon, PA 15642
www.exone.com

ExOne GmbH
Daimlerstrasse 22
86368 Gersthofen
Germany

ExOne KK
161-5 Haneo
Odawara-shi, Kanagawa
Japan 256-0804